
Dental 3D Printer

User Manual

Dear customer,



Thank you for choosing RayForm products!

Maybe you are familiar with 3D printing technology or have purchased RayForm printers before. However, we still highly recommend you read this manual carefully, as the installation techniques and precautions can help you avoid any unnecessary damage or frustration.

1. Please visit <https://www.rayformtech.com> to contact us if you have any questions. You can also learn more information from the website, such as software, videos, models.

2. Social media platforms:

The official social media platform will regularly release updated information on the 3D printing industry, cutting-edge technical news, and common product technical guidelines.

RayForm Official

Copyrighted by "Shenzhen RayForm Technology Co., Ltd ", all rights reserved.
Team RayForm

Safety Instructions

Always follow the safety instructions during assembly and usage, to avoid unnecessary damage to the 3D printer or individual injury.

1. Please contact our Customer Service if you have any issues after receiving the products.
2. In case of emergency, please immediately cut off the power of the 3D printer and contact our technical support.
3. UV light is harmful to the eyes; please avoid direct contact. when operating, wear protective equipment such as anti-UV goggles and gloves.
4. RayForm 3D printer includes components that can cause injury.
5. When the Temperature control resin vat heats up for a long time its surface temperature becomes high. Avoid touching the resin vat without protection after heating.
6. Be cautious when using the scraper, and make sure to orient the sharp parts of machine and tool away from people.
7. Keep the RayForm 3D printer and its accessories out of the reach of children.
8. Vapors or fumes may be irritating at operating temperatures Always use the RayForm 3D printer in an open and well ventilated area.
9. Do not expose RayForm 3D printer to any water or rain environment.
10. Use RayForm 3D printer in an environment with a temperature of 5°C-40°C and a humidity of 20%-80%. For optimal performance, do not exceed this range. Also, avoid direct sunlight exposure.
11. Do not disassemble RayForm 3D printer, please contact technical support if you have any questions.

After-sales service

Process Support: support@rayformtech.com



Customer Service: contact@rayformtech.com


Contents

- 1, Product Overview4
- 2, In the Box5
- 3, Technical Specification6
- 4, Slice File7
- 5, Recommended Print Parameters8
- 6, Start Printing9
- 7, Recycled Resin11
- 8, Machine Maintenance.....12
- 9, FAQ14


Product Overview

In the Box

RF-3024D

Print platform

Power Supply


Manual

USB drive

Metal shovel

Funnel

Protective mask

Protective equipment

Plastic shovel

Technical Specification

Operating System

System	RayForm 3D Printer System
Operation Screen	5-inch capacitive touch screen
Software	ChiTuBox (Support .ctb file)
Connectivity	USB Drive, WLAN

Specifications

LCD Screen	6.6 inch 4K
Light Source	COB light
XY Resolution	4098*2560
Z axis Accuracy	0.01 mm
Suggested Layer Thickness	0.01~0.15 mm
Power Supply	24V - 3A

Physical Dimensions


Dimension	300(L)*280(W)*430(H) mm
Build Volume	A: 143(L)*90(W)*80(H) mm B: 85(L)*85(W)*80(H) mm
Weight	13kg

WLAN

Interface	RJ45 100M
WIFI	2.4Ghz


Slice File

- 1. You can choose to install the slicing software provided in the USB drive or download it from the official website (www.chitubox.com) and install it on your computer.
- 2. Since the RF-3024D is an innovative three-platform integrated 3D printer, when using CHITUBOX Basic, additional device information and print platform parameter settings are required. (Two different platform parameter settings)

Platform	Name	Resolution		Size		
		X	Y	X	Y	Z
	Normal	4098	2560	143	90	80
	Small			85	85	80
	Dual			143	90	80

- 3. Import the 3D model into CHITUBOX Basic, and perform operations such as repair, shell removal, hole drilling and adding supports. The software installation package and operation guide can be found at: <https://docs.chitubox.com/en-US/chitubox-basic/latest/ui-and-features/configure-print-parameters>
- 4. Setting the printing parameters in the software, you can modify the required parameter groups in the resin options. USB drive also have three types of resin models slice files: orthodontic tooth molds, elastic bases, and temporary crowns. Select the appropriate file for printing based on your actual needs.

- 5. After completing the settings, perform the slicing operation and save the sliced file in the .ctb format. Then, import it into the USB drive.


Recommended Print Parameters

Parameters Group	Dental Model Resin	Invisible Denture Base Resin	TEMP C&B Resin
Layer Thickness	0.05 mm	0.05 mm	0.05
Exposure Time	3.5 s	6.5 s	3.7 s
Rest Time After Retract	1.5 s	0 s	0 s
Bottom Exposure Time	20 s	30 s	35 s
Bottom Layers	3	4	4
Lifting Distance	7 mm	8 mm	7 mm
Retract Distance	7 mm	8 mm	7 mm
Lifting Speed	120 mm/min	60 mm/min	90 mm/min
Anti-alias	2		
Notice	1. RF-3024D is suitable for many dental cases. The above data are for reference only, and the specific effect is subject to the actual environment. 2, All RayForm production resins can adapt to the use of equipment, more printing parameters and Settings can visit the official website or consult us. 3. The ambient temperature of resin printing is recommended to be 20°C-30°C(68°F-86°F). If the ambient temperature is too low, the temperature controller can be opened to adjust the temperature of the printing room. Set the PV temperature higher than the SV temperature, greater than 2°C, and the heater will automatically start working. 4. For the first use, please test according to our parameters, do not modify at will. If you need to improve the printing speed, please contact us first. 5. The use of ACF release film can print 45000 layers normally. If you need to replace it, please use the same type.		




--The data above root in RayForm lab, only for reference.

Start Printing

1. Place the equipment on a stable work table, in a well-ventilated environment, avoid strong ultraviolet light, such as fluorescent UV disinfection lamp, so as not to affect the printing effect of abnormal curing.
2. According to the items to be printed, choose the appropriate printing platform and the corresponding material box to install, turn on the power supply, hit the boot key, and pour the required resin.
(Note: for long standing colored resin, shake for 2-3 minutes)

3. Select the appropriate language in the system interface Settings, and then select the pattern for exposure detection.

4. Import the prepared slice file into the U disk, and then insert the U disk into the USB interface on the right of the printer


Suggestions for using U disk:
① It is recommended to use the USB flash drive that comes with the machine. If you use other USB sticks, make sure the USB stick has a capacity of less than 4G and that the format supports FAT/FAT32.
② Please put the print file in the root directory of the U disk to avoid abnormal file reading.

Start Printing


5. Select the model file to be printed on the printer touch screen and enter the detail page.

• The U disk file can be copied to internal storage, and the next printing can be initiated directly on the device.

• Select the files that need to be printed, and then click "Confirm Print".




- **Print completed and lifted to the top:** After the printing is finished, the printing platform automatically rises to the top of the Z-axis.

6. Before printing, printer will conduct an intelligent monitoring to ensure that the printing process can proceed normally. If any error situation are detected during the monitoring process, the machine will automatically report an error, please handle it according to the tips.

Recycled Resin

1. After the printing is completed, wait until the remaining resin on the platform no longer drips. Then remove the platform. Next, use a metal spatula to remove the model and rinse the residual liquid resin on the model surface with 95% alcohol (or other cleaning agents). After cleaning and drying, perform other post-processing steps such as post-curing on the model.

2. When there is some solidified resin in the material box, plastic shovels can be used to clean the material box. Failure to clean in time will affect the subsequent printing process and may even cause damage to the release film and the exposure screen.

3. Remove the material box, pour out the remaining resin from the inclined opening of the box, and filter the resin using a filter screen and funnel. The funnel model file is provided in the complimentary USB drive. If needed, you can print and use it.


4. If the machine is not used for a long time, thoroughly clean the material box and store the resin in a light-proof and sealed container to prevent any deterioration that could affect its next use.

5. If the stored resin has not been used for a long time, it is necessary to shake it for 2-3 minutes before the next use to prevent the precipitation of pigments and other organic substances, which may affect the printing results.


Machine Maintenance

Material box maintenance

- Remove the resin remaining on the release film: Use the tray cleaning function. For the final step, use a plastic spatula to scrape off the solidified resin layer by layer; do not use sharp objects to scrape the release film as it may cause damage.

- Replace the release film: Users can view the cumulative number of printing layers on the device's self-check interface and replace the release film in a timely manner according to the actual situation to avoid affecting the printing quality or causing resin leakage.

Reset the number of print layers after replacement

- If not used within 48 hours, please filter the resin and store it in a light-proof and sealed container.
- If resin gets on the tray, please wipe it clean immediately.

Z-axis maintenance

If there is a frictional noise during the Z-axis operation, please apply an appropriate amount of lubricating grease to the Z-axis screw rod.

Machine Maintenance

Firmware upgrade

- **Local upgrade:** Download the latest firmware from the official website and copy it to the USB drive. Be sure to delete any previous versions of the firmware on the USB drive. Then insert the USB drive into the printer and complete the upgrade.

Machine cleaning

- **Clean the printing platform:** Wipe it clean directly with a tissue or clean it with alcohol.
- **Protect the hardened surface:** If there is any resin adhering to the scratch-resistant film, please replace the scratch-resistant film immediately.
- **Clean the body:** Wipe it clean with alcohol.

FAQ

Complete the initial problem investigation by referring to the following solutions, or contact the after-sales support for additional assistance.

Printing

Model does not stick to the platform

- The exposure time at the bottom is insufficient. Please increase the exposure time.
- The contact area between the model's bottom and the platform is small. A bottom valve needs to be added.

2. Model fracture of the layer

- The machine shakes during the printing process.
- The release film becomes loose after long-term use and needs to be replaced.
- The printing platform or the material box is not tightened.
- The lifting speed is too fast.
- The shell model has not been drilled.

3. The model has layering patterns and deformations

- Check if there is insufficient support.
- Reduce the lifting speed.

4. There are tufts of something resembling kelp attached to the tray or the model.

- Excessive exposure has occurred. Therefore, the exposure time at the bottom and the normal exposure time need to be reduced.

Network connection

1. WLAN connection failed

- The WLAN name or password is incorrect. Reset the network and try reconnecting.
- The WLAN signal is poor or there is no network. Ensure that the current network environment is normal, then reset the network and reconnect.

2. Cannot add printer


- Check "Server Location". For the Chinese mainland region, please select "China"; for other countries and regions, select "Global".

Processing of electronic components

This sign indicates that this equipment must not be disposed of together with household garbage. When the equipment reaches the end of its lifespan, please handle it in accordance with local environmental protection regulations. It can be recycled at the electrical and electronic waste recycling points. By properly disposing of the equipment, the improper handling of electronic waste can be avoided, which would cause harm to the environment and health. Recycling also helps to protect natural resources. Do not discard this product together with unclassified household garbage.

Once again, we would like to thank you for choosing RayForm products! We offer a maximum one-year warranty for the products (and accessories). If you encounter any problems, please log in to the RayForm official website (<https://www.rayformtech.com/>) to inquire about solutions or contact the corresponding store customer service. Our professional after-sales technical team will be dedicated to serving you.

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