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This document aims to:

- **Define critical operational aspects** including equipment precautions, structural components, operating procedures, I/O port testing, software installation, and Rayform software operations;
- **Guide equipment operators** in daily usage and troubleshooting through standardized references;
- **Standardize the operating procedures for the post-curing oven**, serving as an operational guideline for product post-processing.

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This operation manual applies to the RF4030 Series 3D Printers.

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Temperature: Maintain ambient temperature at 24 - 28℃.

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Humidity: Less than 65% RH.

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Ventilation: The air should be well-circulated.

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Lighting: Use safety yellow light (illuminance < 80uw/cm²).

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Power Supply: 220V-50Hz (±10%) AC, single-phase (rated power: 2.5kW)

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Grounding: The equipment must be properly grounded for safety.

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Power Isolation: Turn OFF the power switch and the main power (circuit breaker) before inspection or maintenance.

Warning: If the inspection and maintenance work is in progress and the main power (circuit breaker) is not turned OFF, the machine may start accidentally. In this case, the machine may be damaged and personal safety may be at risk.

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Emergency Response: If smoke, abnormal heating, or unusual noise is detected:

Immediately turn OFF the power switch and the main power (circuit breaker)
Contact our technical support team.

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Replacement Parts: Use only manufacturer-recommended components.

Note: Non-approved parts void the warranty.

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Storage: Store the resin in a designated area, protected away from light and heat, and keep it locked.

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Wear medical-grade gloves during operation.

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Skin Contact: Rinse immediately with soap and water.

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Eye Contact: Flush with clean water or saline solution for 15 minutes.
Seek immediate medical attention.

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Ingestion: Rinse mouth thoroughly and induce vomiting. Seek immediate medical attention.

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Anhydrous Alcohol: For post-print cleaning (excludes water-washable or self-cleaning resins).

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Medical Gloves: Prevent direct contact with resin/alcohol.

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Absorbent papers: Wipe residual resin/alcohol from models.

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Angled Spatula: Safely remove cured models.

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The system comprises five core modules: Structural Components, Mechanical Motion System, Optical Engine Assembly, Hardware Control Unit, and Software Control System.

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The structural assembly provides foundational support for both the optical engine assembly and mechanical motion system.

• It comprises:

1. Equipment Frame: Main chassis (typically aluminum alloy or steel) for rigidity.
2. Support Brackets: Reinforced mounts for motion system integration.
3. Resin Vat: for resin containment.
4. Optical Engine Mounting Platform: Precision-machined surface for DLP module installation.

• Installation Layout:

1. Optical Engine Assembly: Mounted on the Optical Engine Mounting Platform.
2. Mechanical Motion System: Anchored to Support Brackets and Equipment Frame through vibration-damped couplings.

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The mechanical motion system governs the movement of critical functional parts. It comprises:

1. Motors: Stepper motors with closed-loop control
2. Motor Drivers: High-precision drivers for microstepping operation
3. Air Pump: Pneumatic system
4. Build Platform
5. Float Tank: Resin level adjustment system
6. Scraper & Squeegee Blade: Titanium alloy tools for resin leveling

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The optical engine assembly is responsible for photopolymerization imaging in the 3D printing process. It comprises the following components:

1. Four Optical Engine Modules: DLP based UV light sources (wavelength: 405nm $\pm$ 5nm).

2. HDMI Data Cables connect to Computer GPU to transmits sliced image data from the computer to the engines.
3. USB Control Cables connect to Industrial USB Hub to controls exposure timing/light intensity via software commands
4. Optical Engine Power Cables connect to Engine Power Switch
Provides stable power supply

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
The hardware control system comprises three core elements:

1. Chassis Server (Industrial Control Computer – IPC)
Function: Hosts and executes the printer control software.
Data Interface: Manages optical imaging data transmission via HDMI/USB protocols.
2. Display Monitor
Role: Human-Machine Interface (HMI) for operational visualization.
3. Liquid Level Sensor
Function: Provides real-time resin level feedback.

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

The software control system is primarily comprised of Rayform® Software, which serves as the centralized control system for both the mechanical motion components and the optical exposure assembly.


- 1、设备总开关 Main Power Switch
- 2、光机电源 Optical Engine Power Supply
- 3、显示器电源 Display Power Supply
- 4、继电器 Relay
- 5、端子排 Terminal Block
- 6、地线排 Grounding Terminal Block
- 7、DLP 光学引擎 DLP Optical Engine
- 8、工控机 Industrial Control Computer (IPC)
- 9、运动控制卡 Motion Control Card
- 10、光机电源适配器 Opto-Mechanical Power Adapter
- 11、滤波器 Power Filter
- 12、电源 Power Supply Unit (PSU)
- 13、GPIO 控制板 GPIO Control Board
- 14、电机驱动器 Motor Driver

- 1、机械运动电源开关 Mechanical Motion Power Switch
- 2、光机电源开关 Optical Engine Power Switch
- 3、电脑电源开关 Computer Power Switch
- 4、设备复位按钮 System Reset Button
- 5、急停按钮 Emergency Stop Button (ESB)

- 1、设备框架 Frame Assembly
- 2、光学引擎组件 Optical Engine Assembly
- 3、平台升降电机 Platform(Z-Axis) Lift Motor
- 4、平台拖筋组件 Platform Lead Screw Assembly
- 5、铲料组件 Model Ejection Kit
- 6、刮刀 Resin Leveling Blade/Scraper
- 7、打印平台 Build Platform
- 8、导料组件 Material Guiding Component

- 1、补液筒 Resin Refill Canister
- 2、补液蠕动泵 Resin Refill Peristaltic Pump
- 3、浮筒 Float Tank
- 4、液位探测器 Liquid Level Sensor
- 5、树脂温度计 Resin Thermometer
- 6、浮筒电机 Float Tank Level Adjustment Motor

4.3. 设备操作 Equipment Operation

按照配套清单清点配件，并将设备移至设备安装区，设备安全区应达到上述室内条件的要求。然后进行设备安装调试，调试完成后，开启设备，启动软件，进行打印。设备操作包括设备开机及操作步骤、模型打印完取件及后处理、设备关机三个部分 Verify all accessories against the packing list, relocate the equipment to the installation area (ensuring compliance with Section 4.1.1 environmental requirements), complete installation and commissioning, then power on the system and initiate printing. Operations include startup procedures, post-printing model handling, and shutdown protocols.

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

Settings Panel contains four menu bars: File, Settings, Control, and Help.

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1) “File” Menu contains two buttons: Open File, and Exit.


- a. Click “Open File” to select the model files for printing.

- b. Click the “Exit” and the software will close.

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The “Settings” menu bar includes six buttons: “System Configuration”, “Motion Configuration”, “Engine Configuration”, “Image Settings”, “Process Templates”, and “Language Settings”.

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a. Click the "System Configuration" button, and the system configuration dialog box will pop up.

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Click "System Configuration" button

Motion Axis: Maximum 4 axes (factory preset)

Optical Engines: Up to 4 DLP modules (factory preset)

Sensors: Limit 14 units (factory preset)

Initial Print Layer is set according to the printing situation. Generally, printing starts from the first layer of the image.

Print start position is the position where the printing platform begins to print the first layer. (Z-axis initial height)

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
Click “Liquid Level Settings” to set the liquid height level and the leveling step size for liquid level adjustment.

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Dynamic exposure does not require setting.

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
Log verbosity is categorized into three levels (0, 1, 2) and typically does not require user configuration under normal operating conditions. Default setting: Level 1 (recommended for daily use)

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
The build platform ejection position is configured for automatic part removal after printing.

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All the above settings are factory-configured.

To apply user-modified parameters:

1. Click [Apply] after adjusting any configuration.
2. Exit the system configuration dialog box and click the [System Initialize] button.
3. The new parameters will be activated in the next print job.

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Click the “Motion Settings” button to open the Motion Configuration dialog box. This interface controls motion parameters for the printer’s motors.

Device ID: Contains 4 motor channels.

Factory Presets: All parameters are factory-optimized and locked.

Modification Protocol: Contact authorized service personnel for adjustments.

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Click the “Optical Engine Configuration” button to open the configuration dialog. This interface sets parameters for optical engines, which are pre-configured at the factory.

• Critical Note:

• Simulation Mode: Check this box to disable UV light output from the selected engine (for safe testing purposes).

• Modification Protocol: Contact authorized service personnel for parameter adjustments.

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Click the “Image Settings” button to open the image configuration dialog box.

1. Edge Enhancement: Sharpens projected contours to optimize the printing effect..


2. Contour Smoothing (Anti-aliasing): Reduces pixelation for smoother edge transitions.


3.Print File Type:

- Direct STL Printing: Raw model printing without pre-slicing.
- Sliced File Printing: Standard mode using pre-sliced layers.

4.Layer Thickness: Defines Z-axis resolution.

5.Distortion Correction: Compensates for optical lens distortion (pre-calibrated per engine).

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

Click the “Process Template” button to open the process configuration dialog.

Parameter Definitions:

- Step is Exposure Stage which divide printing into multiple phases with distinct parameters. If a single set of process parameters can be applied to print a full dental model, only one Step is required in this configuration
- Image Layer is total number of sliced image layers in the current stage
- exp (s) is exposure time which is UV curing duration per layer (critical for resin polymerization)
- AXISZD is Z-Axis Retract Distance (mm)
- DDELAY is Post-Print Settling Delay (s) which is the time for resin to stabilize after printing
- AXISZU is Z-Axis Post-Retract Lift (mm) which is reserved for advanced dynamic leveling.
- FDELAY is Pre-print Settling Delay (s). It is Resin stabilization time before layer exposure


f. 点击“语言设置”按钮，弹出语言设置对话框。可选择“中文”、“英文”2种语言。
Click the “Language Settings” button to open the language configuration dialog box, where you can select between ‘Chinese’ and ‘English’.

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The “Control” menu includes four functional buttons: “Platform Control”, “Optical Engine Control”, “Functional Test”, and “Vision Calibration”.

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Click the “Platform Control” button to open the motion control dialog box, which includes four modules: Build Platform Control, Squeegee Blade Motion Control, Float Tank Level Adjustment and Scraper Control.

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Click the “Optical Engine Control” button to open the engine control dialog box, which includes two modules: Device Information and Optical Engine Control.

Key Functions:

Individual Exposure Testing:


Engine ID: Select the optical engine (1-4) for standalone light output verification.

Exposure Current: Set the LED driver current .

Click “Current Test” to validate light output.

Prerequisites:

Ensure DLP projection and LED power are enabled before testing.

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
Click the “Functional Test” button to open the diagnostics dialog box, which includes three modules: System Check, Component Test and I/O Test.

I/O Test (Critical Focus)


Verifies the input/output signals of motion components via four test options

1. Brake Control is IOBit1. Toggle the switch to test motor brake engagement.
2. Air Pump Control is IOBit2. Activate/deactivate the pneumatic pump.
3. Liquid supplement control is IOBit3. Open/close the resin refill solenoid valve.
4. Liquid temp control is IOBit4. Enable/disable resin tank temperature regulation.

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The Liquid Level Settings interface utilizes data from the liquid level sensor to achieve continuous real-time calibration of the resin height level.

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Click the “Liquid Level Calibration” button to open the calibration dialog box. This function must be operated by authorized service personnel only.

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
The “Help” menu includes two buttons: “Device Registration” and “About”

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When the device approaches its license expiration date, click the “Help Documentation” button to open the Device Registration dialog box. Follow these steps:

- 1.Generate Authorization Serial Number and save the license key to your desktop.
- 2.Email the key to Shenzhen Rayform company and wait the authorization file.
- 3.Save the received authorization file to your desktop.
- 4.Return to this interface and click “Import Authorization File” to complete activation.

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b. Other Software Interface Elements

4.5.1. Image Display Area is Preview Window

Function: Projects real-time layer images during printing.

4.5.2. Temperature Display is CPU Temperature Monitor

Function: Displays the temperature of the industrial control computer's CPU.

4.5.3. Print Progress Bar

Function: Visualizes current print progress (percentage + layer count).

4.5.4. Time Status Panel

Elapsed Time: Total duration since printing started.

Remaining Time: Estimated time to completion based on current progress.

4.5.5. Resin Level Indicator

Function: Displays real-time resin height level.

4.5.6. Print Control Buttons

a. Start Print icon: Initiates the printing process.

b. Pause/Resume icon: Temporarily halts printing; resume by clicking .

c. Abort Print icon: Terminates the print job irreversibly.

4.5.7. Action Log

Function: Records real-time device operations (e.g., motor movements, exposure events).

4.5.8. Model Metadata Panel

Print Model Details: Display detailed information of the model.

4.5.9. File Selector

Function: Load STL/BMP model files (must enable file type in Settings → Image Settings).

4.5.10. Reconnect/Refresh Button

Function:

1. Re-establishes communication with optical engines and motors.
2. Synchronizes and updates configuration parameters.

4.5.11. Pre-Print Setup Button

Function: Prepares the printer for printing by:

1. Positioning the build platform at the initial Z-height.
2. Adjusting the float tank to set the resin level to the calibrated position.

4.5.12. Manual Motion Control Button

Function: Enables manual adjustment of Build Platform, Squeegee Blade, and Float Tank

4.5.13. Login/Logout Button

Function: Authenticates user access. Logs in and out the current user session.